



www.participant.life
[/participantap](https://www.facebook.com/participantap)
[@participantap](https://twitter.com/participantap)
[@participant_ap](https://www.instagram.com/participant_ap)



**Durable &
Affordable**



Versatile & Reliable

Designed by Participant in California 



Why Pony Lite?

The Pony Lite stands out for durability and quality among dedicated transport wheelchairs.

Research shows that most wheelchairs have their first component failure within six months. Whether that is a broken wheelchair footplate in an elder care facility or a loose wheel lock at a caregiver's home - the wellbeing of the wheelchair rider is diminished.



Versatile

Pony Lite has solid, robust wheels that allow mobility where thinner and lighter wheels may bog down or break.

Frail elders need safety and Pony Lite offers key features:

- Swing away and removable footrests so that tripping risk is minimized.
- Secure and easily adjustable wheel locks.
- Flip back and removable armrests so that transfers to the toilet or bed are unobstructed.



mm (or as noted)

Seat Width	500
Seat Depth	430
Back Height	450
Front Seat Height / Rear Seat Height	465
Seat Angle from the Ground	3 degrees
Footrest Height	280 - 410
Wheelbase	430
Overall Width	680
Overall Length	1,080
Product Weight / Net Weight	19 kg
Weight Capacity	130 kg
Push Handle Height	950
Rear Wheel Dimensions	16 x 1" (400 x 25 mm) wheel with solid PU tire on plastic mag hub
Rear Wheel Camber Angle	0 degrees
Caster Dimensions	200 mm diameter and 50 mm wide solid PU tire
Packaging	One unit per carton; partially assembled; 10 minutes assembly time required.
Package Size	720 x 390 x 660
Package Weight / Gross Weight	22.5 kg
Certifications and Testing	ISO 9001, ISO 13485, CE, Tested ISO 7176 complete, 2X on double drum durability test - 400,000 cycles, ISWP Caster Test - Pass, Conveyor belt track test calibrated to real world conditions - Pass. (see website for all tests and test reports)
Warranty Details	Frame: 3 years Components: 6 months.